


## Experiment No. 7

**Aim:** To carry out the assay of Isonicotinic acid hydrazide (Isoniazid ) Tablets IP.



### Requirements:

### Apparatus:

Volumetric flask, Measuring cylinder, Analytical balance, Weight box, Beaker, Burette, Conical flask and UV Spectrometer.

### Chemicals:

Isoniazid tablets, Potassium bromate, Methyl red, Dimethyl formamide, Sodium carbonate, Sodium thiosulfate, Starch solution.

### Principle:

The reaction is between isoniazid and potassium bromine, but solution of bromine is not stable. Therefore, a small amount of potassium bromide is added to the acidified solution of isoniazid, which is slowly titrated with potassium bromate. The reaction between  $\text{KBrO}_3$ ,  $\text{KBr}$ ,  $\text{HCl}$ , liberate bromine and this bromine oxidizes isoniazid. At the end point when  $\text{HCl}$  get deflected, changes colour from red to yellow.

### Preparation and Standardization of standard solutions:

#### 1. Preparation of Potassium Bromate, 0.0167 M:

Dissolve 2.783 gm of potassium bromate in sufficient water to produce 1000 ml.

#### 2. Preparation of Sodium Thiosulphate, xM:

Dissolve  $248x$  g of Sodium thiosulphate and  $2x$  g of sodium carbonate in sufficient water to produce 1000 ml.

#### 3. Standardization of Potassium Bromate:

Transfer the volume of about 30 ml of potassium bromate into a glass stoppered flask, add 3 ml of potassium iodide, following by 3 ml of hydrochloric acid. Allow to stand for five minutes at room temperature, then titrate the liberated iodine with standard sodium thiosulphate by using 3 ml of starch solution as an indicator. Collect for a blank on the same quantities of the same reagents. Each ml of 0.1 N sodium thiosulphate is equivalent to 0.002784 gm of potassium bromate.

**Procedure:****1. Assay Method by Volumetric Analysis:**

- Weigh and powder 20 tablets.
- Weigh accurately quantity of the powder equivalent to 0.4 g of isoniazid and dissolve in water, filter and wash the residue with sufficient water to produce 250 ml.
- Add 50 ml of water to 50 ml of the resulting solution, 20 ml of hydrochloric acid and 0.2 g of potassium bromide.
- Titrate slowly with 0.0167 M potassium bromate with continuous shaking using 0.05 ml of methyl red solution as indicator, until the red colour disappears.
- Each ml of 0.0167 M potassium bromate is equivalent to 0.003429 g of  $C_6H_7N_3O$ .

**2. Assay Method by UV-Spectroscopic Method:**

- Weigh about 80 mg of isoniazid and add 150 ml of dimethyl formamide, then add sufficient water to produce 500 ml.
- Dilute 5 ml to 100 ml with water and mix. Measure the absorbance of the resulting solution at the maximum at about 367 nm.
- Calculate the content of  $C_6H_7N_3O$  taking 750 as the specific absorbance at 367 nm.

**Result:**

The given sample contains ..... mg of isoniazid.